Spectrochemical Supplies

INCLUDING

ELECTRODES
PLATES
and
FILM

INDEX	PAGE
Serial Marker	1
Glow-Box	1
Spacing Dividers	1
Spectroscopic Plates and Film	2
Photographic Chemicals	2
Plastic Electrode Stand	3
Electrode Funnel and Tweezers	3
Parr Pellet Press	3
Graphite Electrodes, Rods and Powders	4-6
Plate Storage Cabinet	7



INDUSTRIES, INC.

P.O. BOX 798
METUCHEN, N. J. 08840

(201)-549-7144

SPECIALISTS IN SPECTROSCOPY


SERIAL MARKER

FOR NUMBERING NEGATIVES IN:

MICROSCOPY, X-RAY DIFFRACTION, ELECTRON DIFFRACTION, METALLOGRAPHY, NUCLEAR PHYSICS, SPECTROSCOPY AND PHOTOGRAPHY

The Serial Marker photographically assigns a 5-digit number to any light sensitive emulsion. In the darkroom, before development, the film or plate—of any size—is held against a Lucite window through which a number is projected by pressing a switch. The number is advanced by pressing a second button. Resetting, if necessary, is accomplished by rotating the side knob.





Model 3702 has a built in attenuator for controlling the exposure, and is particularly useful in those laboratories where several emulsions, differing in speed, are used. The instrument is first calibrated for each emulsion. This is done by determining the proper position of the attenuator for the emulsion with, say, a 2-second exposure. Once calibrated, the SERIAL MARKER will consistently print a correctly exposed number on similar emulsions. It is equipped with an adjustable gate so that the identifying number may be printed in the same relative place on each film or plate.

3702 Serial Marker, for numbering photographs consecutively. Operates on 115 volts ac; attenuator for adjusting light intensity; size of 5-digit number may be varied from 3 x 15 mm to 4 x 20 mm; fence for aligning plates or films; spare 7.5 watt lamp.

Each \$83.00

GLOW-BOX

Here's a viewing box that should be evaluated in terms of its versatility and usefulness rather than its low cost. The 11" x 9" viewing surface, of 1/4" thick Plexiglas, is illuminated with a high-intensity, rapid-start, cool, circular fluorescent lamp. The overall dimensions of the box (11-2/4" x 11-3/4" x 4" high) permit it to be stored and used in a desk drawer where it is instantly available yet out of the way when not wanted. Retractable legs permit the box to be flush-mounted or tilted for better viewing. The dimensions of the viewing surface are sufficient for examining two 4" x 10" plates or an 8 1/2" x 11" film. Other suggested uses for the Glow-Box are: 1) illuminating samples as they are being titrated; 2) color comparison; 3) tracing; 4) comparing graphs on different sheets of paper.

3710 Glow-Box, viewing box as described above; weight 7 lb, complete with 22-watt, circular fluorescent lamp, 115 vac, carrying handle; finished in baked hammertone grey enamel. 12" x 12" x 4", 11 lb Each \$44.00





3506 Spacing Divider, 6" model, an accurately constructed pantograph-like drafting tool, has eleven teeth so arranged that they always divide the extreme setting into ten equal parts. Each tooth is numbered, left to right on one side, right to left on the other. Constructed of stainless steel, its maximum spread is 23 mm, minimum distance between teeth is 3 mm. This instrument is recommended for most spectrographic wavelength measurements.

Each \$ 43.00

3512 Spacing Divider, 12" model, same as above, but maximum spread 47 mm and minimum distance between teeth 6 mm.

Each \$ 55.00

RECOMMENDED SPECTROGRAPHIC EMULSIONS

SPECTROSCOPIC PLATES AND FILM

EMULSION Spectrum Analysis #1

MOST USEFUL SPECTRAL RANGE

2500-4400A

APPLICATIONS

Moderately slow, fine grained, high contrast and resolving power. Wavelength range minimizes overlapping orders. Most popular emulsion for quantitative and semiquantitative work.

DELIVERED PRICES

	1	FILM	1	PLATES	
SA#1	Daylight SA/413	Darkroom SA/421	Per doz.	Per 3 dz.	Per Case
4" x 10" 2" x 10" Individual Rolls	12.90	12.15	8.40	22.70 14.60	172.50 208.00
6 rolls or more	11.22	10.57	1		
Spectrum Analysis #3	2500-5000A	g	Speed about twice that of SA#1 grain, medium contrast and resolvi concentration range coverage is a	ng power. Excelle	nt when wide
SA#3 4" x 10"	SP/413		9.00	25.30	100.00
Individual Rolls 6 rolls or more	14.00 12.18		5.00	25.30	192.30
#33	2500-4400A	а	This emulsion strikes an average be and is superior to SA#3 in grain background interferes using SA#3.	iness. For trace a	SA#3 in speed analysis when
4" x 10" 2" x 10"	I		1	27.00 18.00	194.00 243.00
103-0**	2500-5000A		Very fast (about 4 times that of SA ditions), moderately coarse grain, power emulsion. Most useful for weak lines or short duration expos	medium contrast analysis of micro	low resolving
103-F**	2500-6750A	I	Extended wavelength coverage for arly useful for Na and Li determina	r 103-0 type emu ations.	Ision. Particu-
1-N**	2500-8750A	i i	Very fast particularly in the reg medium contrast and resolving po analysis of K, Rb and Cs. Despite t should not be used as a general and low speed in the u.v. portion deep-freeze storage is recommende	wer. Particularly its broad wavele emulsion because of the spectrum.	useful for the ngth coverage of granularity
4" x 10" 2" x 10"	SP/417		21.50	_	464.40
Individual Rolls 6 Rolls or more	23.30 20.25		14.20	_	573.75
Q-2 Ilford Plates**	below 2500A	5	Medium speed, grain, contrast and sensitive to handling because of i determination of elements whose	its thin emulsion.	Useful in the


**Since Eastman Kodak recommends keeping these emulsions at temperatures no higher than 55°F, shipments are normally made via the fastest means and are packaged with dry-ice. The cost of premium transportation plus a \$5.00 packaging charge is additional in such cases.

30.40

PHOTOGRAPHIC CHEMICALS

			5 gal.
D-19 Developer	one gallon	\$1.30	
Kodak Fixer	one gallon	.85	
Kodak Rapid Fixer	one gallon	1.75	6.20
Indicator Stop Bath	16 oz. btl.	1.30	





PLASTIC ELECTRODE STAND

3051	Plastic Electrode Stand, for 1/4" electrodes	Each \$18.00
3052	Plastic Electrode Stand, for 3/16" electrodes	Each \$18.00
3053	Plastic Electrode Stand, for 1/8" electrodes	Each \$18.00

ELECTRODE FUNNEL and TWEEZERS

3001 Funnel for filling 1/4" electrodes
3002 Funnel for filling 3/16" electrodes
3003 Funnel for filling 1/8" electrodes

3503 Tweezers, stainless steel, for handling spectroscopic graphite electrodes Each \$ 3.30 Dozen \$31.60

Pair \$ 3.30





PARR PELLET PRESS

(SPEX SPEAKER, Vol VIII-No. 1)

This hand operated pellet press provides enough force against a small area to prepare, from powder, a suitable tablet for the cup of an electrode. The plunger and die set of 0.118" diameter produces pellets of the proper size for L-3909, L-4000 and L-4006 Union Carbide preforms and likewise pellets from the 0.178" diameter die set fit L-3900, L-4012 and L-4018 preforms.

3625 Pellet Press, hand operated, without plunger and die, 9" x 5" x 10", 18 lbs. \$ 93.00

3626 Plunger and Die, for preparing pellets 0.118" dia., fill height of die 3/4" \$ 57.00

3627 Plunger and Die, for preparing pellets 0.178" dia., fill height of die 3/4" \$ 57.00

GRAPHITE ELECTRODES, RODS and POWDERS

Spex has been preeminent as a supplier of graphite products for many years because:

- .. we maintain a large inventory from the world's leading manufacturers
- ...deliveries are made on a same-day basis, fob destination within the continental USA.
- .. "special" preforms are ordered in advance of scheduled shipments and stocked in sufficient quantities to provide immediate deliveries
- .. we "hand-pick" electrodes free of particular contaminants for your special applications
- ..our staff includes experienced spectrochemists who can offer sound advice when it is required.

In order to continue superior service at current prices in the face of rising costs, we are introducing SPEX HIGH PURITY GRAPHITE PRODUCTS.

In doing so we pledge to provide you with the finest graphite products available today in terms of purity, electrical characteristics and physical dimension tolerances, as set by ASTM. Further, we pledge to keep prices, wherever possible, at levels established some years ago. One way to hold the price line is to prepare large quantities of the most popular designs and eliminate those which have shown little, if any, acceptance.

The accompanying list of preformed electrodes, rods and powders is the result of a culling process that took into account sales during the last few years. Many shapes which have been only marginally utilized are not present in this list; we have almost cut in half the variations in our 1966 catalog and still retain most of the ASTM standard preformed electrodes. These are more than ample for the majority of spectrochemical techniques. The **impurity content** of all items is guaranteed to be a **maximum of 6 ppm** for the total of all elements.

WHEN ORDERING, PLEASE NOTE:

- Use Spex catalog numbers for the Spex-labeled graphite products. If your order references only our catalog number you will receive Spex electrodes in a high purity, normal density (HPND) grade. The same catalog numbers with a "D" suffix will bring you the high density (HPHD) grade of high purity electrodes.
- For any electrode or graphite product not listed here, please contact us. We shall happily quote on special orders and quantities.
- Should you prefer the National brand, specify their catalog numbers for each item; we maintain a complete stock of National products for your convenience.



SPECTROSCOPIC PREFORMED ELECTRODES

Cat. N	HPHD#			100 Pcs.		ASTM		No.
						No.	AGKSP+	SPK
			SA	MPLE ELECTRO	DES			
				Necked Crater				
1000	4000D	1/4"		\$22.00	5/32" deep	S-12	L-3912	L-371
042	4042D			22.00	3/16" deep		L-4018	L-421
001	4001D	3/16"	· =	22.00	3/16" deep	S-13	L-3903	L-370
004	4004D		3	19.00	3/32" deep		L-3906	L-370
005	4005D			19.00	3/16" deep	S-14	L-3909	L-370
029	4029D			19.00	3/16" deep		L-4000	L-420
030	4030D		2000000	19.00	3/32" deep		L-4006	L-420
033	4033D	1/8"	-	19.00	1/8" deep		L-3905	L-370
				Crater				
002	4002D	1/4"		20.00	3/16" deep	S-8	L-3900	L-370
016	4016D			20.00	1/16" deep	S-5	L-3982	L-378
020	4020D	1/8"	_	16.00	1/4" deep		L-3979	L-377
034	4034D			16.00	0.059" deep		L-3975	L-377
035	4035D		-	16.00	3/16" deep		L-3977	L-377
				Micro-cup				
1079	4079D	1/4"	-	22.00	0.054"x0.062"			L-425
				Solution Analysi	s			
011*	4011D+	1/2"	0 -	14.00	1/8" thick disc	D-1	L-4075	L-427
1027	4027D		0 =	16.00	0.200" thick disc	D-2	L-4072	L-427
1074	4074D	1/4"		35.00	Vacuum Cup, 3/8"x2"			L-379
1075	4075D			35.00	Vacuum Cup, 5/8"x2"			L-378
074A		3/4"		12.00/ 10	Teflon Cup for 4074 & 4075			L-379
		7120		ode Caps & Ped		202		
1018	4018D	1/4"	-	16.00	9/32" deep, thinwall	S-3	L-3918	L-371
1017	4017D	1/8"		16.00	1-1/2" long pedestal	S-1	L-3919	L-37
Packaged	200/box							

	LABEL	Dia.	Price/	DESCRIPTION		"National"	
Cat IPND+	. No.	- 1	100 Pcs.		ASTM	Cat.	No.
TPNDT	HPHD#				No.	AGKSP+	SPK
			COUNTER ELECTR	ODES			
			Pointed				
1010	4010D	1/4"	\$20.00	120°	C-2	L-3966	L-376
1019	4019D	1/8"	16.00		C-1	L-4036	L-423
		91	Rounded				
1041	4041D	1/4"	20.00	1/16″r	C-5	L-3957	L-375
1009	4009D	3/16"	24.00	1/16"r, double ended, 2" long		L-3955	L-375
1039	4039D	-	17.00	1/16″r		L-3951	L-375
1040	4040D	-	18.00	1/16"r, 2" long		L-3954	L-375
1073	4073D	1/8"	17.00	1/32″r		SP1009	
			Flat				
1024	4024D	1/4"	18.00		C-3	L-3921	L-372
1072	4072D	3/16"	14.50		C-9	L-3923	L-372
1071	4071D	1/8"	14.00		C-6	L-3922	L-372
			GRAPHITE RODS, 12"	LONG**			
005	10050	E (10 × / 00E)			"National"		
085 086	4085D	5/16" (.305) 1/4" (.242)	45.00/25		\$ 92.50/50	L-3812	L-383
087	4086D	3/16" (.180)	35.50/25		71.00/50	L-3809	L-3829
7.5.31	4087D		30.00/25		120.00/100	L-3806	L-382
088	4088D	1/8" (.120)	49.00/50		98.00/100	L-3803	L-3823
'Also avai	lable in 4" a	nd 6" lengths upon					
***	•		ULAR GRADE GRAPHITE F	ODS, 12" LONG			
409		1/4" (.242) 3/16" (.180)	20.00/100		24.00/100	L-4	309
409 409		1/8" (.120)	18.00/100		23.00/100	L-4:	306
409	2	(1.20)	19.00/100		22.00/100	L-43	303
			GRAPHITE POWE	ER			
406	1			-100 mesh, for	Buffering		SP-2>
406				-200 mesh, for	Buffering		SP-2
406	3			-200 mesh, for	Briquetting		SP-10

⁺HPND and AGKSP are 1.61g/cc density graphite #HPHD and SPK are 1.90g/cc density graphite



PLATE STORAGE CABINET

By commercial standards, the ordinary 4" x 10" spectrographic plates are an odd size and do not fit office file cabinets. The pictured unit is ideal for the purpose. It contains 9 drawers, each 11-1/4" x 4-5/8" x 12" deep and equipped with two separators per drawer. About 1000 plates in envelopes can be stored in the unit. The separators as well as the front of the drawers may be labeled for reference.



3820 Plate Storage Cabinet, with nine drawers for 4" x 10" plates. Overall dimensions 19-3/4" high x 34" wide by 12" deep. Steel, welded frame construction, finished in grey baked enamel. 74 lbs.

Each \$ 66.00

ORDERING INFORMATION

Terms are net 30 days to rated firms. To avoid delays purchasers who have not transacted any previous business with Spex Industries should include commercial references or remittances with the initial order.

Shipments will be made F.O.B. Destination anywhere within the continental United States via our choice of surface transportation. Any premium transportation, packaging or invoicing requested will be additionally invoiced.

Minimum order \$10.00



INDUSTRIES, INC. - BOX 798, METUCHEN, N.J. 08840

