

FAX: 201-603-9647

Initial Concentration of SPEX Standard, ppm • 201-549-7144

DILUT-ULATOR  
Final Concentration of Diluted Standard, ppm



STANDARDS

Triple Checked For Quality



DIRECTIONS

1. Locate the Concentration (in ppm) of the SPEX Standard you want to dilute on the "Initial Concentration of Standard Scale A".
2. Align the "A" Concentration with the final Concentration (ppm) to which you want to dilute on the "Final Concentration Scale B".
3. Find the final (total) volume of the Solution you want to make (the size of the volumetric you are using) on the "Final Volume Scale C".
4. Read the Aliquot of original standard required from the "Initial Volume of Standard Scale D".

© SPEX Industries, Inc. 1990




STANDARDS AND COMPOUNDS

3880 Park Avenue, Edison, N.J. 08820

Ph: 201-549-7144 Outside NJ: 800-LAB-SPEX

FAX: 201-603-9647

Triple Checked For Quality



STANDARDS



Final Volume of Diluted Standard, ml  
Part No. D-901

DILUT-ULATOR  
Initial Volume of SPEX Standard Required, ml

Final Volume of Diluted Standard, ml  
Part No. D-901

DILUT-ULATOR  
Initial Volume of SPEX Standard Required, ml

FAX: 201-603-9647

Initial Volume of SPEX Standard Required, ml

AD

BC

CB

DA